

Types of Vaporizers



1. Tec 3 Isoflurane Key



2. Tec 3 Isoflurane Funnel



3. Tec 3 Ex Sevoflurane Funnel



4. Tec 3 Ex Isoflurane Key



5. Aladin Cassette Sevoflurane Key



6. Aladin Cassette Sevoflurane Easy-Fil



7. Tec 4 Isoflurane Funnel



8. Tec 4 Isoflurane Key



9. Tec 4 Isoflurane Bottom Mount



10. Tec 4 Halothane Funnel



11. Tec 4 Ethrane Funnel



12. Ohio Isoflurane Funnel



13. Bickford Isoflurane Funnel



14. Tec 5 Sevoflurane Key



15. Tec 5 Isoflurane Funnel



16. Tec 5 Isoflurane Key



17. Tec 5 Halothane Key



18. Tec 7 Sevoflurane Easy



19. Tec 7 Isoflurane Easy



20. Draeger 19.1 Isoflurane Funnel



21. Draeger 19.1 Sevoflurane Key



22. Draeger 19.3 Isoflurane Key



23. Draeger 2000 Sevoflurane Key



24. Draeger 2000 Isoflurane Key



25. Draeger 2000 Sevoflurane Funnel, Quick or Draeger Fil



26. Draeger 2000 Isoflurane Funnel, Quick or Draeger Fil



27. Universal Drawover



28. Tec 6 Desflurane




29. Penlon Sigma Delta Sevoflurane Funnel




30. Penlon Sigma Elite Sevoflurane Key

EASY FIL

USED FOR: TEC 7



FUNNEL FILLERS/ANTI-SPILL (VET USE ONLY)




DRAEGER FILLS

USED FOR: DRAEGER 2000




KEY FILLERS

USED FOR: TEC 3, TEC 4, TEC 5, OHIO, DRAEGER 19.1, 19.3, 2000, PENLON SIGMA DELTA



QUICK FIL

USED FOR: ABBOTT
DRAEGER 2000



Metal Draeger Key Filler

Any questions call 1-800-717-5955



1.

Dr. 2000 Non-Auto Exclusion (Draeger Comp. Block)




2.

Dr. 2000 Auto Exclusion (Draeger Comp Block)



3.

Dr. 2000 Fixed (Cage Mount)



4.

Dr. 2000 Ohmeda Mount (Ohmeda Mount Comp Block)



5.

Tec 3 Cage Mount



6.

Dr. 19.3-Draeger Back (Draeger Comp. Block)



7.

Dr. 19.3 Ohmeda Back (Ohmeda Comp. Block)



8.

Dr. 19.1 Fixed (Cage Mount)



9.

Tec 4 Bottom Mount



10.

Tec 7 Ohmeda Back, -Looks like back of Tec 4, & 5 (Ohmeda Mount Compatibility Block)



11.

Tec 5 Ohmeda Back with Ohmeda Mount Compatibility Block and Pole Mount


Comp Block

Pole Mount

Fillers on vaporizers



Dräger 19.1 Sevoflurane Key



Key Fill




Sold by GAS

Key Fill Adaptor



Dräger 19.1 Isoflurane funnel fill



Funnel Fill



Sold by GAS

Funnel Fill Adaptor

Dräger 2000's can have 4 different fillers.



Dräger 2000 Sevoflurane



Quick Fil
(Abbott Only)




Quick Fil Adaptor



Draeger Fil



Draeger Fil Adaptor



Funnel Fill
(No Adaptor Needed)



Key Fill



Key Fill Adaptor

Tec 7's can have 2 different fillers.



Tec 7 Isoflurane



Funnel Fill
(No Adaptor Needed)



Easy Fil



Easy Fil Adaptor




Compatibility Block **Ohmeda Mount**



Draeger 2000
Ohmeda Mt.



Tec 4,
Tec 5, Tec 7



Compatibility Block
Pole Mount

Tec 5 Ohmeda
Back with Ohmeda
Mount Compatibility
Block and Pole
Mount



Compatibility Block **Draeger Mount**



Draeger
19.3




Draeger 2000
Auto Exclusion



Draeger 2000
Non-Auto
Exclusion



Compatibility Block **Cage Mount**



Draeger
2000 Fixed



Draeger
19.1 Fixed

Any questions call 1-800-717-5955

Other Parts



Ohmeda Mount
Compatibility Block
with Pole Mount



Pole Mount



23 mm Connector Sets



Silicone Tubing